



Malaria Curative Kit

Anti-malaria medicine 2

This kit is designed to diagnose and rapidly treat an attack of malaria to avoid the risk of death. It should always be kept within reach when one has been exposed to malaria. It will be especially useful when returning from a malarial zone, when a traveler is in a country where malaria is unknown.



DIAGNOSTIC

Contents

- 1 forehead thermometer
- 3 plastic lancets
- 3 white packets
- 3 rapid diagnostic tests for malaria falciparum:
 - 3 immunoquick strips
 - 3 tubes of reagent
- 3 disinfectant wipes

Use of kit

1. Hold a lancet by its sides
2. Turn the movable blue cap to free it
3. Carefully disinfect the tip of a finger
4. Apply the exposed needle to the finger
5. Press trigger (a painless prick)
6. Press finger to produce a drop of blood
7. Apply the blotting test end of strip to the drop of blood
8. Insert strip into tube
9. After 5-15 minutes read the result

TREATMENT

Contents

- 1 box of a complete treatment for an attack of malaria (12 tablets).

Use of medicine

1. Immediately take 4 tablets of this medicine with food or milk.
2. Take 4 more tablets 24 hours later
3. Take 4 tablets at the same time on day 3.

Once you have started the treatment, it must be completed, even if you begin to feel better.



Wherever possible, use your MALARIA CURATIVE KIT in consultation with a doctor. If you are unable to see one rapidly (within 6 hours following the first sign of symptoms), you can perform the malaria diagnostic test yourself, preferably with the assistance of a person accompanying you.



Solutions For Your Good Health